

Date: 13/04/2015

ADAL 3D-Walking analysis


The new tools of biomechanics

Dynamometric treadmill for the continuous dynamic measurement on each leg, and recording of the 3 spatial components (3D) of the ground reaction forces while walking.

Applications

- Normal and pathological walking study (after effects of hemiplegia, neuromotive troubles)
- Assistance in pre-operation diagnosis and post-operation follow up.
- Quantification of functional impotence, measurement of progress.
- Rehabilitation (recovery of functional locomotive independence)
- Dynamic qualification of external prostheses for amputees.
- Clinical research from quantitative criterion.


Measurement equipment for standard model








- Equipped with 2 Kistler 3 axis sensors on each side, left and right. Totally 4 sensors and 6 charge amplifiers.
- 6 analog channels (3 for each side) for the measurement of the following raw data's:
 - Vertical forces (Z axis): Fz right and left range 3 000N
 - Antero posterior forces (Yaxis): Fy right and left range +/- 500N
 - Medio lateral forces (X axis): Fx right and left range +/- 500N
- 2 quick numerical inputs for the instantaneous measurement of the displacement and speed of the belts left and right (from 2 encoders' 3600pts/200mm).
- This configuration allows the measurement of global forces: Fztotal, Fytotal, Fxtotal.
- Delivered with a National Instrument data acquisition board (Ref. PCI 6024E) and the adapted acquisition and analysis software: ADiSOFT2004©.







Date: 13/04/2015

Main features:

Belt speed 0 to 10km/h (6.2 miles)

Walking speed 0.1 to 10km/h (0.3 to 6.2 miles)

Motor drives 2 brushless motors with feedback encoders

Main power 3Kw (2*1.5kW)

Supply voltage 200 to 240Vac /16A / 50-60Hz – Single phase with earth grounding

Emergency stop Pushbutton on control box and front barrier

Progressive stop on the control box

Adjustable side and front hand rails.

Speed control:

- Manually with push button on the control box,

- From the ADISOFT software (separately on both side if needed)


- From other devices by RS232 connection

Software: ADISOFT2000©

This software's, easy to use, thanks to its base modules (monitoring, acquisition, analysis and report) enables an objective and detailed gait analysis.

Running under Windows® on any PC type computer, enables:

- Measurement and evaluation of tridirectional ground reaction forces for each leg separately
- Vertical, antero posterior, medio lateral force analysis with automatic calculation of its characteristic points
- Analysis of step series (min. and max average, average deviation).
- Calculation of symmetry index (Herzog and Kaufman methods)
- A report can be printed; data can be exported under Excel ®....
- Language French, English, other languages on request.





Date: 13/04/2015

Other possibilities:

Can be installed on the floor or in a pit, in this case the belt and the floor are at the same level.

Installed on the floor, Belts altitude 290mm Here installed and synchronized with a LOCOMAT robot.

The Locomat robot is a product of The Switzerland Company HOCOMA www.hocoma.com





Installed in a pit, Belts and floor at the same level Here connected to a motion analysis system VICON.

www.vicon.com



The equipment is delivered with a calibration kit including a strain gauge sensor, a mobile electronic and adapted pieces to introduce known forces in horizontal and vertical axis.

